

Product datasheet for GA211940

Mouse Wdr90 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr90 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 90
Symbol:	Wdr90
Components:	GA211940G1, Wdr90 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211940G2, Wdr90 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211940G3, Wdr90 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZPG2</u>
Synonyms:	3230401M21Rik; AI551153; mKIAA1924
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)