

Product datasheet for **GA211943**

Mouse Vmac activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Vmac CRISPRa kit - CRISPR gene activation of mouse vimentin-type intermediate filament associated coiled-coil protein
Symbol:	Vmac
Components:	GA211943G1, Vmac gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211943G2, Vmac gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211943G3, Vmac gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BP01</u>
Synonyms:	AI662250
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)