


Product datasheet for **GA211944**

Mouse Cyp4f15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp4f15 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 4, subfamily f, polypeptide 15
Symbol:	Cyp4f15
Components:	GA211944G1, Cyp4f15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211944G2, Cyp4f15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211944G3, Cyp4f15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	AI787289; Cyp4f-15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)