


## Product datasheet for **GA211948**

### Mouse Catsperd activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | Catsperd CRISPRa kit - CRISPR gene activation of mouse cation channel sperm associated auxiliary subunit delta                                                                                                                                                                                                      |
| Symbol:       | Catsperd                                                                                                                                                                                                                                                                                                            |
| Components:   | GA211948G1, Catsperd gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211948G2, Catsperd gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211948G3, Catsperd gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>E9Q9F6</u>                                                                                                                                                                                                                                                                                                       |
| Synonyms:     | 619718; 4921529N20Rik; 4933402B14Rik; AW045609; Gm6095; Tmem146                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)