

Product datasheet for GA211961

Mouse Arap3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arap3 CRISPRa kit - CRISPR gene activation of mouse ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3
Symbol:	Arap3
UniProt ID:	<u>Q8R5G7</u>
Synonyms:	Centd3; cnt-d3; Drag1; E030006K04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)