

Product datasheet for **GA211967**

Mouse Psmg2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psmg2 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) assembly chaperone 2
Symbol:	Psmg2
Components:	GA211967G1, Psmg2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211967G2, Psmg2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211967G3, Psmg2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9EST4</u>
Synonyms:	1700017I17Rik; AW545363; Clast3; Tnfsf5ip1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)