

## Product datasheet for GA211987

### Mouse Tm9sf3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tm9sf3 CRISPRa kit - CRISPR gene activation of mouse transmembrane 9 superfamily member 3
Symbol:	Tm9sf3
UniProt ID:	<u>Q9ET30</u>
Synonyms:	1810073M23Rik; 2810031D16Rik; AI115521; AI413748; AW146116; AW549777; mKIAA4036; Smbp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)