

## Product datasheet for **GA212022**

### Mouse Pi4kb activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pi4kb CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol 4-kinase beta
Symbol:	Pi4kb
UniProt ID:	<u>Q8BKC8</u>
Synonyms:	ESTM41; PI4K-beta; PI4Kbeta; Pik4cb
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)