

Product datasheet for **GA212072**

Mouse **Chrn3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Chrn3 CRISPRa kit - CRISPR gene activation of mouse cholinergic receptor, nicotinic, beta polypeptide 3
Symbol:	Chrn3
UniProt ID:	<u>Q8BMN3</u>
Synonyms:	5730417K16Rik; Acrb3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)