

Product datasheet for GA212115

Mouse Slc35b3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc35b3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 35, member B3
Symbol:	Slc35b3
UniProt ID:	Q922Q5
Synonyms:	4921526O06Rik; AI428480; CGI-19; PABST2; PAPST2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

