

## Product datasheet for **GA212120**

### Mouse Rnf187 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Rnf187 CRISPRa kit - CRISPR gene activation of mouse ring finger protein 187                                                                                                                                                                                                                                  |
| Symbol:       | Rnf187                                                                                                                                                                                                                                                                                                        |
| Components:   | GA212120G1, Rnf187 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212120G2, Rnf187 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212120G3, Rnf187 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8BFX1</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | 2410016F01Rik; AB030190; RACO-1                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)