

Product datasheet for **GA212126**

Mouse Cops8 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Cops8 CRISPRa kit - CRISPR gene activation of mouse COP9 signalosome subunit 8 |
| Symbol: | Cops8 |
| UniProt ID: | <u>Q8VBV7</u> |
| Synonyms: | 9430009J09Rik; AA408242; Csn8; Sgn8 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

