


Product datasheet for **GA212139**

Mouse Pnrc1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pnrc1 CRISPRa kit - CRISPR gene activation of mouse proline-rich nuclear receptor coactivator 1
Symbol:	Pnrc1
Components:	GA212139G1, Pnrc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212139G2, Pnrc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212139G3, Pnrc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	5730463C12Rik; B4-2; Prol2; PRR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)