


## Product datasheet for **GA212146**

### Mouse Flacc1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Flacc1 CRISPRa kit - CRISPR gene activation of mouse flagellum associated containing coiled-coil domains 1                                                                                                                                                                                                    |
| Symbol:       | Flacc1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA212146G1, Flacc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212146G2, Flacc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212146G3, Flacc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8BVM7</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | 4933405P16Rik; 4933425F06Rik                                                                                                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)