

Product datasheet for **GA212180**

Mouse Pfdn4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pfdn4 CRISPRa kit - CRISPR gene activation of mouse prefoldin 4
Symbol:	Pfdn4
Components:	GA212180G1, Pfdn4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212180G2, Pfdn4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212180G3, Pfdn4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	0610025M10Rik; 1810063B06Rik; C-1; C1; Pfd4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)