

Product datasheet for **GA212210**

Mouse Tspan9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspan9 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 9
Symbol:	Tspan9
Components:	GA212210G1, Tspan9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212210G2, Tspan9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212210G3, Tspan9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BJU2</u>
Synonyms:	6720474K14Rik; 9430079M16Rik; AU018597
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)