

Product datasheet for **GA212233**

Mouse Tbx10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tbx10 CRISPRa kit - CRISPR gene activation of mouse T-box 10
Symbol:	Tbx10
Components:	GA212233G1, Tbx10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212233G2, Tbx10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212233G3, Tbx10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q810F8</u>
Synonyms:	Dc; Tbx7; Tbx13
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)