

Product datasheet for GA212245

Mouse Hyal3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hyal3 CRISPRa kit - CRISPR gene activation of mouse hyaluronoglucosaminidase 3
Symbol:	Hyal3
UniProt ID:	Q8VEI3
Synonyms:	Hyl3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)