

Product datasheet for **GA212275**

Mouse Syt17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Syt17 CRISPRa kit - CRISPR gene activation of mouse synaptotagmin XVII
Symbol:	Syt17
Components:	GA212275G1, Syt17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212275G2, Syt17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212275G3, Syt17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q920M7</u>
Synonyms:	Bk
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)