

Product datasheet for **GA212291**

Mouse **Slc35b1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Slc35b1 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 35, member B1
Symbol:	Slc35b1
Components:	GA212291G1, Slc35b1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212291G2, Slc35b1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212291G3, Slc35b1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P97858</u>
Synonyms:	Ugalt2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)