

## Product datasheet for **GA212361**

### Mouse **Kcnq3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	Kcnq3 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily Q, member 3
Symbol:	Kcnq3
UniProt ID:	<a href="#">Q8K3F6</a>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)