

Product datasheet for **GA212364**

Mouse *Scn4a* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scn4a CRISPRa kit - CRISPR gene activation of mouse sodium channel, voltage-gated, type IV, alpha
Symbol:	Scn4a
UniProt ID:	Q9ER60
Synonyms:	mH2; Nav1.4; SkM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)