

Product datasheet for GA212395

Mouse Vmn1r43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Vmn1r43 CRISPRa kit – CRISPR gene activation of mouse vomeronasal 1 receptor 43
Symbol:	Vmn1r43
Components:	GA212395G1, Vmn1r43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212395G2, Vmn1r43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212395G3, Vmn1r43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VIC9</u>
Synonyms:	V1ra2; V1ra5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

