

## Product datasheet for GA212403

### Mouse Vmn1r25 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | Vmn1r25 CRISPRa kit – CRISPR gene activation of mouse vomeronasal 1 receptor 25                                                                                                                                                                                                                                  |
| <b>Symbol:</b>       | Vmn1r25                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>   | GA212403G1, Vmn1r25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212403G2, Vmn1r25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212403G3, Vmn1r25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | Vlrc8                                                                                                                                                                                                                                                                                                            |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

