

Product datasheet for **GA212423**

Mouse Tle6 activation kit by CRISPRa

Product data:

| | |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Tle6 CRISPRa kit - CRISPR gene activation of mouse transducin-like enhancer of split 6 |
| Symbol: | Tle6 |
| Components: | GA212423G1, Tle6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212423G2, Tle6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212423G3, Tle6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | Q9WVB3 |
| Synonyms: | 1810057E06Rik; Grg6 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)