

Product datasheet for **GA212428**

Mouse Ly6g5c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ly6g5c CRISPRa kit - CRISPR gene activation of mouse lymphocyte antigen 6 complex, locus G5C
Symbol:	Ly6g5c
Components:	GA212428G1, Ly6g5c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212428G2, Ly6g5c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212428G3, Ly6g5c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K1T5</u>
Synonyms:	G5c; NG33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)