

## Product datasheet for **GA212471**

### Mouse Olfr74 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Olfr74 CRISPRa kit – CRISPR gene activation of mouse olfactory receptor 74                                                                                                                                                                                                                                    |
| <b>Symbol:</b>       | Olfr74                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA212471G1, Olfr74 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212471G2, Olfr74 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212471G3, Olfr74 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | mOR-EV; MOR174-4                                                                                                                                                                                                                                                                                              |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

