

## Product datasheet for **GA212476**

### Mouse Pip4k2c activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pip4k2c CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-5-phosphate 4-kinase, type II, gamma
Symbol:	Pip4k2c
UniProt ID:	<a href="#">Q91XU3</a>
Synonyms:	Pip5k2c
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)