

Product datasheet for **GA212478**

Mouse Ttyh2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttyh2 CRISPRa kit - CRISPR gene activation of mouse tweety family member 2
Symbol:	Ttyh2
Components:	GA212478G1, Ttyh2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212478G2, Ttyh2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212478G3, Ttyh2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TH73</u>
Synonyms:	1110001A03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)