

Product datasheet for GA212483

Mouse Stk33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Stk33 CRISPRa kit - CRISPR gene activation of mouse serine/threonine kinase 33
Symbol:	Stk33
UniProt ID:	Q924X7
Synonyms:	4921505G21 Rik; AW061264
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)