

Product datasheet for **GA212512**

Mouse **Kcnn3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Kcnn3 CRISPRa kit - CRISPR gene activation of mouse potassium intermediate/small conductance calcium-activated channel, subfamily N, member 3
Symbol:	Kcnn3
UniProt ID:	P58391
Synonyms:	KCa2.3; SK3; SKCA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)