

Product datasheet for GA212517

Mouse Eri3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Eri3 CRISPRa kit - CRISPR gene activation of mouse exoribonuclease 3
Symbol:	Eri3
Components:	GA212517G1, Eri3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212517G2, Eri3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212517G3, Eri3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C460</u>
Synonyms:	PINT1; Prnpip1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)