

Product datasheet for GA212520

Mouse Plxnb3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plxnb3 CRISPRa kit - CRISPR gene activation of mouse plexin B3
Symbol:	Plxnb3
UniProt ID:	Q9QY40
Synonyms:	AI451018; Plxn6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)