

## Product datasheet for **GA212537**

### Mouse P2ry14 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2ry14 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled, 14
Symbol:	P2ry14
UniProt ID:	<a href="#">Q9ESG6</a>
Synonyms:	A330108O13Rik; Gpr105; P2Y14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)