

Product datasheet for GA212570

Mouse Usp48 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Usp48 CRISPRa kit - CRISPR gene activation of mouse ubiquitin specific peptidase 48
Symbol:	Usp48
Components:	GA212570G1, Usp48 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212570G2, Usp48 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212570G3, Usp48 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3V0C5</u>
Synonyms:	2810449C13Rik; AI115503; BC021769; D330022K21Rik; Usp31
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)