

Product datasheet for GA212574

Mouse Papln activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Papln CRISPRa kit - CRISPR gene activation of mouse papilin, proteoglycan-like sulfated glycoprotein
Symbol:	Papln
Components:	GA212574G1, Papln gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA212574G2, Papln gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA212574G3, Papln gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q9EPX2</u>
Synonyms:	E030033C16Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

