

## Product datasheet for **GA212586**

### Mouse *Kcnh7* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh7</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 7
Symbol:	<i>Kcnh7</i>
UniProt ID:	<a href="#">Q9ER47</a>
Synonyms:	9330137111Rik; erg3; Kv11.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)