

Product datasheet for GA212618

Mouse Ms4a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ms4a3 CRISPRa kit - CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 3
Symbol:	Ms4a3
UniProt ID:	Q920C4
Synonyms:	HTm4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)