

Product datasheet for GA212643

Mouse *Asic3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Asic3</i> CRISPRa kit - CRISPR gene activation of mouse acid-sensing (proton-gated) ion channel 3
Symbol:	<i>Asic3</i>
UniProt ID:	Q6X1Y6
Synonyms:	<i>Accn3</i> ; AW742291; DRASIC; SLNAC1; TNAC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)