

## Product datasheet for **GA212650**

### Mouse Vmn1r234 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | Vmn1r234 CRISPRa kit - CRISPR gene activation of mouse vomeronasal 1 receptor 234                                                                                                                                                                                                                                   |
| Symbol:       | Vmn1r234                                                                                                                                                                                                                                                                                                            |
| Components:   | GA212650G1, Vmn1r234 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212650G2, Vmn1r234 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212650G3, Vmn1r234 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | V1rf1                                                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)