


Product datasheet for **GA212765**

Mouse Tet3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tet3 CRISPRa kit - CRISPR gene activation of mouse tet methylcytosine dioxygenase 3 |
| Symbol: | Tet3 |
| Synonyms: | B430006D22Rik; BC037432; D230004J03Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)