

Product datasheet for **GA212787**

Mouse Trim40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim40 CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 40
Symbol:	Trim40
Components:	GA212787G1, Trim40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212787G2, Trim40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212787G3, Trim40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UWA4
Synonyms:	Gm35; Gm319; Gm865
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)