

Product datasheet for **GA212800**

Mouse **Ccdc154** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc154 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 154
Symbol:	Ccdc154
Components:	GA212800G1, Ccdc154 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212800G2, Ccdc154 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212800G3, Ccdc154 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6RUT8</u>
Synonyms:	Ccdc157; Gm317
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)