

Product datasheet for **GA212818**

Mouse Thsd4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Thsd4 CRISPRa kit - CRISPR gene activation of mouse thrombospondin, type I, domain containing 4
Symbol:	Thsd4
Components:	GA212818G1, Thsd4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212818G2, Thsd4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212818G3, Thsd4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UTY6</u>
Synonyms:	ADAMTSL6; Adamtsl4; ADAMTSL6; AI585901; B230114P05Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)