

Product datasheet for **GA212838**

Mouse Vcpkmt activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Vcpkmt CRISPRa kit - CRISPR gene activation of mouse valosin containing protein lysine (K) methyltransferase
Symbol:	Vcpkmt
Components:	GA212838G1, Vcpkmt gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212838G2, Vcpkmt gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212838G3, Vcpkmt gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8C436
Synonyms:	Gm71; Mettl21d; Mettl221d; VCP-KMT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)