

Product datasheet for GA212840

Mouse Pif1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pif1 CRISPRa kit - CRISPR gene activation of mouse PIF1 5'-to-3' DNA helicase
Symbol:	Pif1
UniProt ID:	Q80SX8
Synonyms:	4631410M14; AI449441
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)