

## Product datasheet for **GA212851**

### Mouse *Tmprss7* activation kit by CRISPRa

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	<i>Tmprss7</i> CRISPRa kit - CRISPR gene activation of mouse transmembrane serine protease 7
<b>Symbol:</b>	<i>Tmprss7</i>
<b>Components:</b>	GA212851G1, <i>Tmprss7</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212851G2, <i>Tmprss7</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212851G3, <i>Tmprss7</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<a href="#">Q8BIK6</a>
<b>Synonyms:</b>	B230219I23Rik; Gm1748
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

