

Product datasheet for **GA212857**

Mouse Ankrd33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ankrd33 CRISPRa kit - CRISPR gene activation of mouse ankyrin repeat domain 33
Symbol:	Ankrd33
Components:	GA212857G1, Ankrd33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212857G2, Ankrd33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212857G3, Ankrd33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BXP5</u>
Synonyms:	A930021G21Rik; Panky; Panky-a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)