

## Product datasheet for **GA212864**

### Mouse Iqcj activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Iqcj CRISPRa kit - CRISPR gene activation of mouse IQ motif containing J                                                                                                                                                                                                                                |
| Symbol:       | Iqcj                                                                                                                                                                                                                                                                                                    |
| Components:   | GA212864G1, Iqcj gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212864G2, Iqcj gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA212864G3, Iqcj gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8BPW0</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | D230050A05; EG208426                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)