

Product datasheet for **GA212901**

Mouse Thnsl1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Thnsl1 CRISPRa kit - CRISPR gene activation of mouse threonine synthase-like 1 (bacterial)
Symbol:	Thnsl1
Components:	GA212901G1, Thnsl1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212901G2, Thnsl1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212901G3, Thnsl1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BH55
Synonyms:	AW413632
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)